



Designed by clinicians
for clinicians and their patients

The D R Burton PōcketPEP®

Oscillating Positive Expiratory Pressure

The PocketPEP® is a handheld OPEP device for lung expansion and promotes secretion clearance. PocketPEP® creates oscillating positive expiratory pressure when the patient exhales through the device, helping to loosen and remove bronchial secretions associated with:

- Pulmonary emphysema
- Chronic Obstructive Pulmonary Disease (COPD)
- Atelectasis
- Chronic Bronchitis
- Bronchiectasis
- Cystic Fibrosis
- Asthma
- Nonproductive cough
- Smoker's cough



PōcketPEP®

PōcketPEP® Features and Benefits:

- Ideal for patient use at home
- The palm-sized PocketPEP® is the smallest OPEP available, less than half the size of other products
- In acute settings, OPEP therapy is intended to reduce the risk of hospital acquired pneumonia. When patients continue use of the convenient OPEP therapy at home, it can reduce risk of pulmonary complications post discharge
- With only one moving part, the PocketPEP® is easy for patients to operate and clean.
- Meets Cystic Fibrosis cleaning guidelines
- Optional T-Piece adapts the PocketPEP® for use with most nebulizers

PocketPEP®

*Designed by clinicians
for clinicians and their patients*




PocketPEP® Specifications:

- PocketPEP, item #530, 1 case of 10 each
- T-Piece adapts PocketPEP to most nebulizers, item #550, 1 case of 50 each
- PocketPEP with T-Piece, item #533, 1 case of 10 each
- PVC free
- Not manufactured from latex natural rubber
- Clean with soapy water, dishwasher, boiling or microwave steam bags
- Cold sterilize with isopropyl alcohol or hydrogen peroxide



Rx only. Single patient, multi use.
Do not share or use on multiple patients.



D·R Burton

INSPIRE FOR LIFE

D R Burton Healthcare LLC
3936 South Fields Street
Farmville, NC 27828

Email: inquiries@drburtonhpi.com
Phone: (252) 228-7038