

**Designed by clinicians
for clinicians and their patients**




PocketPEP®

The PocketPEP® is a handheld OPEP device for lung expansion and promotes secretion clearance. PocketPEP® creates oscillating positive expiratory pressure when the patient exhales through the device, helping loosen and remove bronchial secretions associated with:

- Pulmonary Emphysema
- Chronic Obstructive Pulmonary Disease (COPD)
- Atelectasis
- Chronic Bronchitis
- Bronchiectasis
- Cystic Fibrosis
- Asthma
- Nonproductive cough
- Smoker's cough

PocketPEP® Features and Benefits:

- Ideal for patient use at home
- The palm-sized PocketPEP® is the smallest OPEP available, less than half the size of other products
- In acute settings, OPEP therapy is intended to reduce the risk of hospital acquired pneumonia
- When patients continue use of the convenient OPEP therapy at home, it can reduce risk of pulmonary complications post discharge
- PocketPEP® is easy for patients to operate and clean
- Optional T-piece adapts the PocketPEP® to accessories



PocketPEP®

Designed by clinicians for clinicians and their patients

Item	Description	Case Size
530	PocketPEP® OPEP Device	10
533	PocketPEP® OPEP Device with T-piece	10
535	PocketPEP® Advantage OPEP Device with IS Adapter	10

- PVC free
- Not manufactured from latex natural rubber



Rx only. Single patient, multi use.
Do not share or use on multiple patients.



D • R Burton Healthcare LLC
3936 South Fields Street
Farmville, NC 27828

EMAIL: inquiries@drburtonhealthcare.com
PHONE: (252) 228-7038
WEBSITE: www.drburtonhealthcare.com

© 2022 D • R Burton Healthcare LLC. All rights reserved.

PocketPEP 713-0822-03

