




Oscillating Positive Expiratory Pressure



Oscillating Positive Expiratory Pressure

The vPEP[®] HC is a handheld OPEP device for lung expansion and promotes secretion clearance. vPEP[®] HC creates oscillating positive expiratory pressure when the patient exhales through the device, helping to loosen and remove bronchial secretions associated with:

- Pulmonary Emphysema
- Chronic Obstructive Pulmonary Disease (COPD)
- Atelectasis
- Chronic Bronchitis
- Bronchiectasis
- Cystic Fibrosis
- Asthma
- Nonproductive cough
- Smoker's cough



vPEP[®] HC Features and Benefits:

- Compact, ergonomic design
- Easily adjustable resistance
- Easy assembly and cleaning with 4 parts, helping to reduce risk of infection
- Ideal for use in the home or hospital

vPEP® HC Specifications:

- Cleaning-only 4 simple parts
- Drain port
- Adjustable resistance
- PVC free
- Not manufactured from latex natural rubber
- Clean with soapy water

PN	Description	Quantity
512	vPEP	case10
512-T	vPEP with T-piece adapter	case10



Rx only. Single patient, multi use.
Do not share or use on multiple patients.



D • R Burton Healthcare LLC
3936 South Fields Street
Farmville, NC 27828

EMAIL: inquiries@drburtonhealthcare.com
PHONE: (252) 228-7038
WEBSITE: www.drburtonhealthcare.com

© 2022 D • R Burton Healthcare LLC. All rights reserved.

VPEP 714-0822-03

