



The Versatile **OxyPAP**[®]

For Positive Airway Pressure (PAP) therapy




- The OxyPAP device is indicated for the treatment and prevention of atelectasis. It also has the ability to provide supplemental oxygen when used with compressed oxygen.
- Provides therapeutic airway pressures between 5 - 20 cm H₂O.



Versatile components of the OxyPAP may be used to deliver oxygen flow and PAP therapy during OPEP therapy

- The component design of the OxyPAP allows the clinician to tailor therapy to individual patient needs.
- Ideal for atelectatic patients with secretion retention.



Ordering Specifications	Item	Description	Case Size
	562	OxyPAP with Manometer and Tubing	10
	563	OxyPAP with tubing	10
 	562-TS	OxyPAP with Manometer, Tubing, and Mask Adapter	10
	562-6		60
	563-TS	OxyPAP with Tubing and Mask Adapter	10

D • R Burton Healthcare LLC 3936 South Fields Street Farmville, NC 27828
 EMAIL: inquiries@drburtonhealthcare.com PHONE: (252) 228-7038 WEBSITE: www.drburtonhealthcare.com

© 2022 D • R Burton Healthcare LLC. All rights reserved. OxyPAP, OxyJET, vPEP, PocketPEP and D • R Burton Healthcare are registered trademarks of D • R Burton Healthcare

OxyPAP 720-0822-03

